

## RECORD SYMPTOMS

### Follow your physician's instructions carefully.

The Patient Diary is an important part of this procedure. Use this diary to document any symptoms you experience.

Questions? Contact your physician or consult [bardydx.com](http://bardydx.com).

**Traveling through airports:** Inform screening personnel that you are wearing the Carnation Ambulatory Monitor before going through scanner. Bring this Patient Diary to show security personnel.

### Security Screening Statement

This person is wearing a Carnation Ambulatory Monitor (records ECG). It was applied under the direction of a physician. The device is currently adhered to the patient's chest to monitor the heart.



Bardy Diagnostics, Inc.

US Customer Service:  
(844) 422-7393

**EC REP**

Medical Product Service GmbH  
Borngasse 20  
35619 Braunfels  
Germany  
EU Customer Service:  
+31 (0)46 7630422

### Symptoms include:

- |                                |                        |
|--------------------------------|------------------------|
| 1) Palpitations                | 5) Shortness of Breath |
| 2) Dizziness / Lightheadedness | 6) Exercise-Related    |
| 3) Fainted                     | 7) Other               |
| 4) Chest Discomfort / Pain     |                        |

Date/Time:	<input type="text" value="05"/> / <input type="text" value="25"/> / <input type="text" value="19"/>	<input type="text" value="6"/> :	<input type="text" value="37"/>	<input checked="" type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom:	<input type="text" value="2"/>	Other:	<input type="text"/>		

Date/Time:	<input type="text"/> / <input type="text"/> / <input type="text"/>	<input type="text"/> :	<input type="text"/>	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom:	<input type="text"/>	Other:	<input type="text"/>		

Date/Time:	<input type="text"/> / <input type="text"/> / <input type="text"/>	<input type="text"/> :	<input type="text"/>	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom:	<input type="text"/>	Other:	<input type="text"/>		

Date/Time:	<input type="text"/> / <input type="text"/> / <input type="text"/>	<input type="text"/> :	<input type="text"/>	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom:	<input type="text"/>	Other:	<input type="text"/>		

Date/Time:	<input type="text"/> / <input type="text"/> / <input type="text"/>	<input type="text"/> :	<input type="text"/>	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom:	<input type="text"/>	Other:	<input type="text"/>		

Date/Time:	<input type="text"/> / <input type="text"/> / <input type="text"/>	<input type="text"/> :	<input type="text"/>	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom:	<input type="text"/>	Other:	<input type="text"/>		

Date/Time:	<input type="text"/> / <input type="text"/> / <input type="text"/>	<input type="text"/> :	<input type="text"/>	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom:	<input type="text"/>	Other:	<input type="text"/>		

Date/Time:	<input type="text"/> / <input type="text"/> / <input type="text"/>	<input type="text"/> :	<input type="text"/>	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom:	<input type="text"/>	Other:	<input type="text"/>		

Date/Time:	<input type="text"/> / <input type="text"/> / <input type="text"/>	<input type="text"/> :	<input type="text"/>	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom:	<input type="text"/>	Other:	<input type="text"/>		

Date/Time:	<input type="text"/> / <input type="text"/> / <input type="text"/>	<input type="text"/> :	<input type="text"/>	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom:	<input type="text"/>	Other:	<input type="text"/>		

Date/Time:	<input type="text"/> / <input type="text"/> / <input type="text"/>	<input type="text"/> :	<input type="text"/>	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom:	<input type="text"/>	Other:	<input type="text"/>		

Date/Time:	<input type="text"/> / <input type="text"/> / <input type="text"/>	<input type="text"/> :	<input type="text"/>	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom:	<input type="text"/>	Other:	<input type="text"/>		

Date/Time:	<input type="text"/> / <input type="text"/> / <input type="text"/>	<input type="text"/> :	<input type="text"/>	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom:	<input type="text"/>	Other:	<input type="text"/>		

## Patient Diary

PATIENT NAME: \_\_\_\_\_


PHYSICIAN NAME: \_\_\_\_\_

HOSPITAL/CLINIC: \_\_\_\_\_

### DEVICE APPLIED (Required):

Prescribed Wear Time:  Days

Date/Time:  /  /   :   AM  PM



### Wearing the Carnation Ambulatory Monitor (CAM):

1. Following your normal daily activities, wear the CAM for the amount of time prescribed by your physician.

2. If you feel symptoms that may be related to your heart, gently push the button and record date/time in this diary. Do **not** push button repetitively or forcefully.



**CAUTION:** If your CAM becomes dislodged, contact your physician for assistance.

**PUSH ONLY ONCE**  
when you feel  
**SYMPTOMS**

## CLINICAL STAFF USE ONLY:



Register the patient information into BDxCONNECT

Patient Name: \_\_\_\_\_

Patient ID: \_\_\_\_\_

Date of Birth: \_\_\_\_\_ Gender: \_\_\_\_\_


Physician Name: \_\_\_\_\_

Hospital/Clinic: \_\_\_\_\_

### DEVICE APPLIED: (Required)

Prescribed Wear Time:  Days

Date/Time:  /  /   :   AM  PM





## INSTRUCTIONS FOR WEARING THE CARNATION AMBULATORY MONITOR

## RECORD SYMPTOMS

More space for recording symptoms →

**Indications** (Select all that apply)

- Palpitations
- Syncope
- Near-Syncope
- Dizziness / Lightheadedness
- Stroke / TIA
- Atrial Fibrillation (AF) / Atrial Flutter (AFL)
- Bradycardia
- Chest Pain
- Shortness of Breath / Dyspnea
- Arrhythmia Therapy Evaluation
- Other: \_\_\_\_\_

**Underlying Heart Disease** (Select all that apply)


- None Known
- Known AF
- CAD
- CHF
- Non-Ischemic Cardiomyopathy
- Congenital Heart Disease
- Other: \_\_\_\_\_

Pacemaker / ICD      YES   NO

- 3.** Wear the CAM at all times, including showering.
  - ⚠ **CAUTION:** Avoid showering, bathing, or exercise immediately following application, and thereafter avoid activities or environments that result in excessive perspiration, as this may result in a decreased period of monitoring.
  - ⚠ **CAUTION:** Strenuous exercise and activities, such as hot yoga or sauna that may result in excessive sweating, should be avoided.
  - ⚠ **CAUTION:** It is normal for the CAM adhesive material to swell in humid environments or when exposed to moisture. Allow the adhesive to dry following activities such as shower or exercise. If desired, gently pat with a dry towel, but do not attempt to reposition the CAM.
  - ⚠ **CAUTION:** Submersion (such as during swimming or bathing) is not advised. Keep showers brief and the CAM out of the direct stream of water.
  - ⚠ **CAUTION:** Skin irritation. Minor itching or irritation is normal. Only remove before the amount of time prescribed by your physician if irritation from the adhesive is not tolerable. Mark date/time of removal, and make note in this Patient Diary.
  - ⚠ **CAUTION:** Poor contact of the CAM with the skin can negatively affect monitoring performance. Secure the CAM back in place if it becomes loose or detached.

- 5.** Wipe off any remaining adhesive on the skin with the Adhesive Remover Wipe Pad provided.
- 6.** Returning the CAM:  
Place the Carnation Ambulatory Monitor and this Patient Diary inside the box.
- 7.** Seal the mailer with the Mailer Sticker and return per the instructions of your physician/nurse.



**Mailer Sticker (Step 7)**

For Patient Frequently Asked Questions please visit  
[www.bardyd.com](http://www.bardyd.com)

**Symptoms include:**

- 1) Palpitations
- 2) Dizziness / Lightheadedness
- 3) Fainted
- 4) Chest Discomfort / Pain
- 5) Shortness of Breath
- 6) Exercise-Related
- 7) Other

Date/Time: 05 / 25 / 19	61 : 07	<input checked="" type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: 2	Other: _____		
Date/Time: MM / DD / YY	HH : MM	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: #	Other: _____		
Date/Time: MM / DD / YY	HH : MM	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: #	Other: _____		
Date/Time: MM / DD / YY	HH : MM	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: #	Other: _____		
Date/Time: MM / DD / YY	HH : MM	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: #	Other: _____		
Date/Time: MM / DD / YY	HH : MM	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: #	Other: _____		
Date/Time: MM / DD / YY	HH : MM	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: #	Other: _____		

Date/Time: MM / DD / YY	HH : MM	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: #	Other: _____		
Date/Time: MM / DD / YY	HH : MM	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: #	Other: _____		
Date/Time: MM / DD / YY	HH : MM	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: #	Other: _____		
Date/Time: MM / DD / YY	HH : MM	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: #	Other: _____		
Date/Time: MM / DD / YY	HH : MM	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: #	Other: _____		
Date/Time: MM / DD / YY	HH : MM	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: #	Other: _____		
Date/Time: MM / DD / YY	HH : MM	<input type="checkbox"/> AM	<input type="checkbox"/> PM
Symptom: #	Other: _____		