


AI for AF Detection

Generate superior diagnostic yield with a P-wave centric monitor



Published Poster Summary

Atrial signal clarity and rhythm specificity are critical when developing artificial intelligence (AI) to distinguish atrial fibrillation (AF) from rhythms that mimic AF.

Authors: Dan Paris, MSCS; Rodney Boleyn, MSEE; Ezra Dreisbach, BSCS; Chuck Dulken, BSEE; Susan Blancher, NP; Melinda Wang, BSN; and Gust H. Bardy, MD from Bardy Diagnostics, Inc. – Seattle, WA

Our validation was comprised of two-hour excerpts of CAM ECG data chosen at random from 50 AF-positive patients and 50 AF-negative patients (200 hours total). AF was diagnosed if variable P-wave morphology was present for at least 30 seconds. AF presence and duration were confirmed by a team of experienced cardiac electrophysiology clinicians. Disagreements between the three validating electrophysiology clinicians were adjudicated at weekly review meetings.

Results

In the validation cohort, there were 58 patients presenting with normal sinus rhythm (NSR), 27 with atrial flutter (AFL), 24 with ventricular tachycardia (VT), 9 with atrioventricular blocks (AVB), and all patients exhibited AT. The AI differentiates AF not only from NSR but also from myriad other arrhythmias, including those that have R-R interval variability such as AFL and AT. Our results were 96.82% sensitive and 99.86% specific with a positive predictivity of 99.79% for detecting 30 seconds of AF or longer.

Conclusions


Our P-wave centric continuous ECG monitoring technology allows our neural network, or AI, to differentiate between AF and a host of rhythms that mimic AF. AI systems that do not make these distinctions may mislead both patients and clinicians.

The CAM patch AI uses low-amplitude, low-frequency content to distinguish AF not just from normal sinus, but from other arrhythmias that mimic AF.

This includes arrhythmias that have R-R variability like atrial flutter with variable conduction, atrial tachycardia with variable conduction, bigeminy, trigeminy, quadrigeminy, and the various types of heart blocks.



Atrial Flutter with Variable Conduction



Atrial Tachycardia to Atrial Fibrillation

Rhythm specificity matters

The difference between a patient **actually having AF** and having an arrhythmia that **mimics AF** can play out in medical and procedural management and stroke risk.

Gain confidence in AF detection with high signal fidelity and low-amplitude, low-frequency P-wave clarity.

One AI output every half second



CAM patch AI is engineered to give an output for every half second of recording, identifying if a patient has AF and showing where AF starts and stops.

Accuracy vs. validation datasets

For AF episodes of 30 seconds or longer

Sensitivity
96.82%

Specificity
99.86%

Positive predictivity
99.79%



To view the 5-min poster presentation, visit

<https://info.bardyd.com/AFYND>



For more clinical evidence, go to www.bardyd.com